



Features:

roe = recombination cloning, over-expression

1. vermilion as selectable marker
2. attB/P site specific integration
3. Amp resistant
4. For over-expression
 - 10X UAS
 - HSP70 promoter
 - Two insulators
 - one ftz intron
5. To simplify cloning
 - Recombination method
 - Two steps, first step clone PCR fragment into Entry vector (PENTR/D TOPO or mENTRY); second step is recombination with this vector. Detailed procedure for small scale cloning can be found on the TRiP website: <http://flyrnai.org/TRiP-HOME.html>.

PLEASE NOTE: To amplify this vector you will need to use One Shot ccdB Survival 2 T1R Competent Cells (Invitrogen Cat. No. A10460), or other gateway competent cells.

664-2543, vermilion; 2568- 2935, attB; 2946-3377, gypsy; 3394-3428, loxp; 3433-3545, 5XUAS
 3551-3584, loxp; 3591-3700, 5XUAS; 3703-3966, Hsp70; 3974-5715, attR; 5749-5929, ftz intron
 5931-6634, SV40 polyA tail; 6636-7073, gypsy

Accession number: GU931386

Jian-Quan Ni (jni@genetics.med.harvard.edu)

