



Features:

1. vermilion is the selectable marker
2. attB/P site specific integration
3. Amp resistant
4. For high expression of the hairpin
 - 10X UAS
 - HSP70 promoter
 - Two introns: one from white gene, another from ftz gene
5. Detailed procedure for small scale cloning can be found on the TRiP website:
<http://flyrnai.org/TRiP-HOME.html>.

664-2543, vermilion; 2568-2935, attB; 2946-2980, loxp; 2986-3096, 5XUAS; 3102-3135, loxp; 3143-3252, 5XUAS
 3262-3518, Hsp70; 3547-3620, white intron; 3651-3833, ftz intron; 3838-4537, SV40

Accession number: GU931383

Ni et al.,2008. Nature Methods, 5: 49-51.

